



APP-075  
Revision 1  
November 2023

## Nanotrap<sup>®</sup> Microbiome Combined 10 mL Automated Protocol using MagMAX<sup>™</sup> Kit and the KingFisher<sup>™</sup> Apex

**Objective:** This protocol uses Nanotrap Microbiome A Particles and Nanotrap Microbiome B Particles and Nanotrap Enhancement Reagent 3 to capture and concentrate microbes in environmental water samples. It is optimized for microbe capture from 10 mL samples and is compatible with MagMAX Wastewater Ultra Nucleic Acid Isolation Kit.

### Materials and equipment:

Sample Type	
Environmental Water Samples	
Concentration Reagent	Vendor
Nanotrap Microbiome A Particles	Ceres Nanosciences; SKU# 44202
Nanotrap Microbiome B Particles	Ceres Nanosciences; SKU# 65202
Nanotrap Enhancement Reagent 3 (ER3) <sup>1</sup>	Ceres Nanosciences; SKU# 10113
Extraction Kit	Vendor
MagMAX Wastewater Ultra Nucleic Acid Isolation Kit	Thermo Fisher Scientific <sup>™</sup> ; Cat# A52606
Materials/Equipment	Vendor
KingFisher <sup>™</sup> Apex with 96 DW Head	Thermo Fisher Scientific; Cat# 5400930
KingFisher Apex 24 Combi head	Thermo Fisher Scientific; Cat# 24079940
KF Apex 96 KF heating block	Thermo Fisher Scientific; Cat# 24075920
KF Apex 24 DW heating block	Thermo Fisher Scientific; Cat# 24075940
KingFisher 24 Deep-well Plate, Barcoded	Thermo Fisher Scientific; Cat#95040470B
KingFisher 24 Deep-Well Tip Comb & Plate, Barcoded	Thermo Fisher Scientific; Cat#97002610B
KingFisher 96 Deep-well Plate, Barcoded	Thermo Fisher Scientific; Cat# 95040450B
KingFisher 96 Plate (200 µL), Barcoded	Thermo Fisher Scientific; Cat# 97002540B
General Reagents	Vendor
Ethanol	VWR; 1006-012
Molecular grade water	VWR; 45001-044

<sup>1</sup> Precipitate can form in ER3 if stored below room temperature. Allow ER3 to return to room temperature to dissolve the precipitate (can invert 2-3 times to help resuspend, do not heat).

## Capture and Extract Microbes using Nanotrap Microbiome Particles

### Procedure:

#### 1. Nanotrap Combined MagMAX KingFisher Apex Procedure-Part 1

1. *Prepare* "Sample Plates 1" and "Sample Plates 2"
  1. Invert environmental water sample 5 times to mix. After inverting, place on a flat surface for 45 seconds.
  2. Add 4,800  $\mu\text{L}$  of environmental water sample to one well (one well per sample) of a new KingFisher 24 Well Deep Well Plate.
  3. Add another 4,800  $\mu\text{L}$  of environmental water sample to the same well location on a second KingFisher 24 Well Deep Well Plate.
    - a) For example, if you loaded a sample into well A1 of the first plate, load the second volume of that sample into well A1 of the second plate.
  4. Add 50  $\mu\text{L}$  of Nanotrap Enhancement Reagent 3 (ER3) Solution to each sample on the two KingFisher 24 Well Deep Well sample plates (100  $\mu\text{L}$  total).
  5. Add 75  $\mu\text{L}$  of Nanotrap Microbiome A Particles to each sample on the two KingFisher 24 Well Deep Well sample plates (150  $\mu\text{L}$  total).
  6. Add 75  $\mu\text{L}$  of Nanotrap Microbiome B Particles to each sample on the two KingFisher 24 Well Deep Well sample plates (150  $\mu\text{L}$  total).
2. *Prepare* "Lysis Plate"
  1. Add 500  $\mu\text{L}$  of MagMAX Microbiome Lysis Solution to a new (the third) KingFisher 24 Well Deep Well Plate matching the number and location of the "Sample Plate" wells.
3. *Prepare* "Tip Comb Plate"
  1. Insert a new tip comb into a new KingFisher 24 Well Deep Well Plate.
4. *Run NT Script (Request file at sales@ceresnano.com)*
  1. Run **10mL\_NTAB\_MagMAX\_24\_Apex\_56\_v2.kfx**
  2. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.
5. Once the protocol is completed, the "Lysis Plate" will contain lysate that is ready to proceed to Part 2 (**\*caution\* sample may be hot**).

#### 2. Nanotrap Microbiome Combined MagMAX KingFisher Apex Procedure-Part 2

1. *Prepare* "MagMAX Bead Binding" Plate
  1. To a new KingFisher 96 Deep Well Plate, add 400  $\mu\text{L}$  of the cleared lysate (NT lysate) from each well of the lysis plate used in "Part 1 step 5" of the protocol. Keep track of which well contains which sample in this new bead binding plate.

2. Add 530  $\mu$ L of MagMAX Binding Solution to each well in which lysate was added.
  3. Add 10  $\mu$ L of MagMAX Proteinase K to each well in which lysate was added.
  4. Add 20  $\mu$ L of MagMAX DNA/RNA Binding Beads to each well in which lysate was added. The total final volume should be 960  $\mu$ L in each sample-containing well of this plate.
  2. *Prepare "Wash Plate 1"*
    1. Add 1 mL of MagMAX Wash Buffer to a new KingFisher 96 Deep Well Plate matching the number and location of the KingFisher 96 Deep Well Plate- *MagMAX™ Bead Binding Plate* wells.
  3. *Prepare "Wash Plate 2"*
    1. Add 1 mL of 80% EtOH to a new KingFisher 96 Deep Well Plate matching the number and location of the KingFisher 96 Deep Well Plate- *MagMAX™ Bead Binding Plate* wells.
  4. *Prepare "Elution Plate"*
    1. Add 100  $\mu$ L of MagMAX Elution buffer to a new KingFisher 96- 200  $\mu$ L plate matching the number and location of the *MagMAX Bead Binding Plate* wells.
  5. *Prepare "Tip Plate"*
    1. Insert the KingFisher 96 Deep Well Comb into a new KingFisher 96 Deep Well Plate
  6. Run Extraction Script (Request file at sales@ceresnano.com)
    1. Run **MagMAX\_96.kfx**
    2. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.
3. Once the protocol is completed, the KingFisher 96-Elution Plate contains eluates that are ready for downstream analysis or can be stored at  $-80^{\circ}$  C.  
*Note: Multiple freeze-thaw cycles may cause degradation.*

---

**Attachments: 2**

*KingFisher™ Apex*

1. *10mL\_NTAB\_MagMAX\_24\_Apex\_56\_v2.kfx*
2. *MagMAX\_96.kfx*